

Website: www.wmvets.co.uk Email: partners@wmvets.co.uk

Use of Equine Inhaler and care of your Spacer device

These instructions accompany your Equihaler/ Baby Spacer device

- Warm the inhaler to body temperature before use
- Shake the inhaler for 30 seconds before first use
- Remove the cap from the inhaler
- Discard the first puff if hasn't been used for more than 12 hours
- Shake the inhaler for 5 seconds to load the dose
- Attach the inhaler to the spacer, keep vertical!
- Place spacer over nostril ensuring a tight seal
- Watch your horse breathe actuate the inhaler just before/ as the horse breathes in
- Allow a full breath before removing the device
- If possible obstruct the other nostril to enforce breath holding for 5-10 seconds, and then allow the horse to breath normally
- Wait at least 30 seconds before shaking the inhaler again
- Shake for 5s before the next puff
- Repeat for the prescribed number of doses

Wash your spacer device weekly in detergent but do not rinse with water

MACCLESFIELD
 VETERINARY
 HOSPITAL
 38 Cumberland Street
 Macclesfield

Cheshire SK10 1BZ

Tel: 01625 501500

Fax: 01625 612240

CONGLETON BRANCH 18 Moody Street Congleton Cheshire CW12 4AP Tel: 01260 273222 WILMSLOW BRANCH 19 Hawthorn Lane Wilmslow Cheshire SK9 5DD Tel: 01625 524422

☐ FARM CENTRE The Barn - Holly Tree Farm Holmes Chapel Road Lower Withington Macclesfield Cheshire SK11 9DT Tel: 01477 571000 WHALEY BRIDGE
Block B
Ringstones Industrial Estate
Whaley Bridge
High Peak
Derbyshire SK23 7PD
Tel: 01663 732564

EQUINE CENTRE
Somerford Park
Holmes Chapel Road
Somerford
Cheshire CW12 4SW
Tel: 01260 280800